

GC Procedure 2008

Turn on the GC before leaving for the field (follow these steps)

- 1) Start Helium flow – Open only horizontal knob to the right of the regulator – regulator on the left should read near the red mark (80 psi) and two carrier gas dials on the side of the GC should read at marks.
- 2) Make sure the following readings are set correctly:
 - a. Inj./Det. Temp. = 110
 - b. Column Temp. = 45
 - c. Attenuation = 16
- 3) Turn on GC
- 4) Set current to 160 – turn up slowly, wait 1 second between clicks.
- 5) Start recorder

Check the tip and change if needed.

Run Helium –

- 1) Attach one end of long tubing to helium tank, turn on a light helium flow.
- 2) Attach other end of tubing to the GC intake.
- 3) Charge
- 4) Open 1
- 5) Open 2
- 6) Allow helium to flow for twenty seconds.
- 7) Close 1
- 8) Wait three seconds.
- 9) Close 2
- 10) Discharge
- 11) Start (Let sample run at least five minutes).
- 12) Detach tube (Make sure helium is turned off).
- 13) Stop

Check uV reading – should be between -1000 and +5000 uV. If not, use zero knobs on GC to fix.

Zero – press Command, N, then Enter. Press Monit to return to screen and check that it did zero.

Check slope and record it on data sheet - check slope by pressing command, SLP Test, and then Enter. After 20 seconds slope will be reported. Record on data sheet.

Prepare DIC syringes

- 1) Rinse syringe three times with the sample.
- 2) Submerge tip of syringe and slowly draw up about 20mL of the sample.
- 3) Tap out air bubbles and expel the liquid until about 12 mL remains.

☀ **Remember to write on the printout which samples or standards corresponds to each set of peaks.**

Run Gas Standards

- 1) Rinse syringe once with standard CO₂
 - a. To fill syringe attach twist end of syringe to tank. Open green knob and close it. Open white knob and allow gas to fill syringe. Remember to close white knob.
- 2) Fill syringe with standard.
- 3) Attach syringe to GC
- 4) Charge
- 5) Open 1
- 6) Open 2
- 7) Inject sample
- 8) Close 1
- 9) Wait three seconds
- 10) Close 2
- 11) Discharge
- 12) Start
- 13) Repeat five times (you can insert the next sample when the previous sample begins to peak) Make sure dials read Charge, Open, Open then insert sample Close 1 wait three seconds Close 2, and Discharge.
- 14) Stop (when all samples have completely eluted).
- 15) Record values
- 16) Rerun any outliers (divide lowest number by the highest number. Correlation must be within 95%).

Run Air and PCO₂ Samples

- 1) Inject 10 mL of room air to clear out standards and start run (charge, open 1, open 2, inject, close 1, wait 3s, close 2, discharge, start)
- 2) Attach first Air syringe to GC
- 3) Charge
- 4) Open 1
- 5) Open 2
- 6) Inject 10 mL of air sample.
- 7) Close 1
- 8) Wait 3 seconds.
- 9) Close 2
- 10) Discharge
- 11) Start
- 12) Repeat for second syringe.
- 13) Stop
- 14) Set timer for five minutes.
- 15) Rerun with remaining 10 mL if the correlation is less than 95%.
- 16) Follow above procedures for pCO₂ syringes (except step 1)

Run DIC Samples

- 1) Expel water until volume is 10mL
- 2) Inject 200 μL of 2N H_2SO_4 directly into tip of syringe.
- 3) Immediately inject He gas into the syringe until the total syringe volume is 30 mL.
- 4) Cap the syringe.
- 5) Repeat steps 1-3 for 2 more syringes
- 6) Shake the capped syringes vigorously for 1 minute.
- 7) Attach 1st syringe to GC tip
- 8) Charge
- 9) Open 1
- 10) Open 2
- 11) Inject gas – about 10 mL (be careful not to inject water into the GC).
- 12) Close 1
- 13) Wait three seconds
- 14) Close 2
- 15) Discharge
- 16) Start
- 17) Repeat steps 7-15 for other syringes.
- 18) Stop (When all peaks have eluted).
- 19) Wait 5 minutes from last peak before starting next run
- 20) Run additional replicates if correlation is less than 95%.

Rerun Gas Standards (same as above).

Check Slope

Run Helium (same as above).

Turn Off GC:

- 1) Slowly turn Current to zero.
- 2) Turn off recorder
- 3) Turn off GC
- 4) Shut off He flow (make sure both tanks are off).

Write the Date and Name of Lake on GC printout.